


## EXPECT RESULTS

**NeoRevive®** Improvement of pebbled appearance



**NeoRevive®** Improvement of fine lines



**NeoBright®** Whitening of brown spots



**NeoBright®** Texture Improvement and skin brightening



CLINICALLY PROVEN | SAFE AND EFFECTIVE



FOR YOUR PERSONAL CONSULTATION CONTACT:

**BBJ**  
  
 01525 719653  
[www.beautybyjanice.com](http://www.beautybyjanice.com)

In the Media



CLINICALLY PROVEN | SAFE AND EFFECTIVE



EXPECT  
BEAUTIFUL  
SKIN

With  
**TriPollar®**  
TECHNOLOGY

geneO<sup>TM</sup>+

4-IN-1  
SUPER FACIAL

CLINICALLY PROVEN | SAFE AND EFFECTIVE

## Geneo+™ Super Facial

Geneo+™ is an advancement in skin care that delivers 4 essential treatments simultaneously:

- Exfoliation of the outer layer of the skin
- Infusion of unique revitalizing formulas
- Oxygenation generated from within
- Tightening & Contouring

### Geneo+™ '4 IN 1' APPROACH

Geneo+ takes its inspiration from natural hot springs - known for inducing skin oxygenation. By simulating this effect, Geneo+ triggers a body response that sends oxygen to the treated area. At the same time Geneo+ exfoliates the skin and creates an optimal environment for infusion of essential nutrients both during and after treatment. Geneo+ is suitable for all skin types and facial areas.


## RADIANT SKIN IS A TREATMENT AWAY

- Visible results from the first treatment
- Skin rejuvenation and firmed complexion
- Oxygenation generated from within
- Optimal absorption of revitalizing nutrients
- Suitable for all skin types

CLINICALLY PROVEN | SAFE AND EFFECTIVE


## HOW IT WORKS

### Exfoliate



Geneo+ exfoliates the upper skin layer to remove dead cells and prepare the skin to receive active nutrients.

### Infuse



Through our exclusive Capsugen, Geneo+ cleanses the skin while infusing nutrient-rich formulas.

### Oxygenate




Geneo+ produces countless minute CO<sub>2</sub> bubbles which gently burst on the skin surface. The Geneo+ applicator is designed to degas the CO<sub>2</sub> bubbles from the gel and infuse them into the epidermis. The applicator's frequency and amplitudes are specifically designed to:

- Produce a kinetic force that optimises the chemical reaction generating CO<sub>2</sub> bubbles.
- Enable optimal penetration of CO<sub>2</sub> bubbles into the epidermis.

The results in optimal absorption of Pollogen's formulas.

### Tighten & Contour

Radiofrequency is used to stimulate fibroblast activity to produce collagen formulation and immediately tighten collagen fibres for tissue retraction in the skin.



CLINICALLY PROVEN | SAFE AND EFFECTIVE

## SKIN CARE TREATMENTS

Geneo+ Technology was designed to work with specially developed nutrient-rich formulas for optimal results.

### NeoRevive®

NeoRevive is an exclusive formula that delivers superior skin rejuvenation and anti-aging.

### NeoBright®

NeoBright provides optimal skin lightening, nourishment and rejuvenation.

## Frequently Asked Questions

### How soon after treatment can I return to my daily activity?

Geneo+ treatments do not require any downtime so you may resume your daily routine immediately. We do, however, recommend using a hypoallergenic sunscreen for at least 24 hours following your session.

### Is the Geneo+ suitable for all skin types?

Yes, Geneo+ is designed to be used on all skin types.

### How long do Geneo+ treatments take?

The treatment takes approximately 55 minutes including pre and post treatment preparation.

## TESTIMONIALS



"Thank you for my gorgeous geneo+ facial. I'm a NEW woman!!"  
Nicola Stephenson - Actress



"Fantastic results with geneo plus, facial was incredible. Thank you"  
Michelle Collins - Actress

CLINICALLY PROVEN | SAFE AND EFFECTIVE